
















 Round bodied needle – taper point – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR20	20 mm		1.5	5/0	75 cm	brown	N17	CC50HR20
			2	4/0	75 cm	brown		CC40HR20
			3	3/0	75 cm	brown		CC30HR20
			3.5	2/0	75 cm	brown		CC20HR20
HR26	26 mm		2	4/0	75 cm	brown		CC40HR26
			3	3/0	75 cm	brown		CC30HR26
			3.5	2/0	75 cm	brown		CC20HR26
HR30	30 mm		3.5	2/0	75 cm	brown		CC20HR30
HR40	40 mm		3.5	2/0	75 cm	brown	N17	CC20HR40
			4	0	75 cm	brown		CC00HR40
			5	1	75 cm	brown		CC01HR40
			6	2	75 cm	brown		CC02HR40
HR48	48 mm		6	2	75 cm	brown		CC02HR48

 Round bodied needle – taper point – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR30	30 mm		4	0	75 cm	brown		CC00DR30
DR45	45 mm		5	1	75 cm	brown		CC01DR45

 Reverse cutting needle – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		1	6/0	75 cm	brown		CC60DS12-75
			1.5	5/0	75 cm	brown		CC50DS12-75
DS19	19 mm		1.5	5/0	75 cm	brown	N17	CC50DS19
			2	4/0	75 cm	brown	N17	CC40DS19
			3	3/0	75 cm	brown	N17	CC30DS19
DS37	37 mm		4	0	75 cm	brown	N17	CC00DS37