







Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	3.5	0	75 cm	green	N18	PE00WN150
	4	1	75 cm	green	N18	PE01WN150
	5	2	150 cm	green	N18	PE02WN150
	2	3/0	250 cm	green	N18	PE30WN250
	3	2/0	250 cm	green	N18	PE20WN250






Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16 16 mm 	1.5	4/0	75 cm	green	N18	PE40HR16
HR20 20 mm 	2	3/0	75 cm	green	N17	PE30HR20
HR26 26 mm 	2	3/0	75 cm	green	N18	PE30HR26
	3	2/0	75 cm	green	N17	PE20HR26
HR30 30 mm 	3	2/0	75 cm	green	N18	PE20HR30
	3.5	0	75 cm	green		PE00HR30
HR36 36 mm 	3.5	0	75 cm	green		PE00HR36
	4	1	75 cm	green		PE01HR36
HR40 40 mm 	4	1	75 cm	green		PE01HR40
	5	2	75 cm	green		PE02HR40




Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS19 19 mm 	1.5	4/0	75 cm	green		PE40DS19
	2	3/0	75 cm	green		PE30DS19
DS24 24 mm 	2	3/0	75 cm	green		PE30DS24
	3	2/0	75 cm	green		PE20DS24
DS30 30 mm 	3	2/0	75 cm	green		PE20DS30
	3.5	0	75 cm	green		PE00DS30
	4	1	75 cm	green		PE01DS30



Round bodied needle – tapercut point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18 18 mm 	1.5	4/0	75 cm	green	N17	PE40HRT18
	2	3/0	75 cm	green	N17	PE30HRT18