






Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	7	3	150 cm	brown		CC03WN150
	8	4	150 cm	brown		CC04WN150





Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR20 20 mm 	1.5	5/0	75 cm	brown		CC50HR20
	2	4/0	75 cm	brown		CC40HR20
	3	3/0	75 cm	brown		CC30HR20
HR26 26 mm 	2	4/0	75 cm	brown		CC40HR26
	3	3/0	75 cm	brown		CC30HR26
	3.5	2/0	75 cm	brown		CC20HR26
HR30 30 mm 	3.5	2/0	75 cm	brown		CC20HR30
HR40 40 mm 	3.5	2/0	75 cm	brown		CC20HR40
	4	0	75 cm	brown		CC00HR40
	5	1	75 cm	brown		CC01HR40
	5	1	90 cm	brown		CC01HR40-90
	6	2	75 cm	brown		CC02HR40
HR48 48 mm 	6	2	75 cm	brown		CC02HR48








Round bodied needle – taper point – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR30 30 mm 	4	0	75 cm	brown		CC00DR30-75
DR45 45 mm 	5	1	75 cm	brown		CC01DR45-75



Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12 12 mm 	1	6/0	75 cm	brown		CC60DS12-75
	1.5	5/0	75 cm	brown		CC50DS12-75
DS19 19 mm 	1.5	5/0	75 cm	brown		CC50DS19
	2	4/0	75 cm	brown		CC40DS19
	3	3/0	75 cm	brown		CC30DS19
DS26 26 mm 	3	3/0	75 cm	brown		CC30DS26
	3.5	2/0	75 cm	brown		CC20DS26
DS30 30 mm 	3.5	2/0	75 cm	brown		CC20DS30
DS37 37 mm 	4	0	75 cm	brown		CC00DS37

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
2	4/0	100 m	brown		CC40100
2	4/0	50 m	brown		CC4050
3	3/0	100 m	brown		CC30100
3	3/0	50 m	brown		CC3050
3.5	2/0	100 m	brown		CC20100
3.5	2/0	50 m	brown		CC2050
4	0	100 m	brown		CC00100
4	0	75 m	brown		CC0075
4	0	50 m	brown		CC0050
5	1	75 m	brown		CC0175
5	1	50 m	brown		CC0150
6	2	50 m	brown		CC0250
7	3	50 m	brown		CC0350
8	4	30 m	brown		CC0430
8	4	25 m	brown		CC0430
9	5	30 m	brown		CC0530
10	6	30 m	brown		CC0630