










 Round bodied needle – taper point – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16	16 mm		1.5	4/0	75 cm	blue		NY40HR16
HR20	20 mm		2	3/0	75 cm	blue		NY30HR20
HR26	26 mm		3	2/0	75 cm	blue		NY20HR26
HR30	30 mm		2	3/0	75 cm	blue		NY30HR30
			3	2/0	75 cm	blue		NY20HR30
HR40	40 mm		3	2/0	75 cm	blue		NY20HR40
			3.5	0	75 cm	blue		NY00HR40
			4	1	75 cm	blue		NY01HR40
			5	2	75 cm	blue		NY02HR40

 Round bodied needle – taper point – 1/2 circle – loop								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR40	40 mm		3.5	0	150 cm	blue	loop	NY00HR40-L
HR50	50 mm		4	1	150 cm	blue	loop	NY01HR50-L

 Reverse cutting needle – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.5	7/0	75 cm	black		NY70DS12
			0.7	6/0	75 cm	black		NY60DS12
			1	5/0	75 cm	black		NY50DS12
			1.5	4/0	75 cm	blue		NY40DS12
DS16	16 mm		0.7	6/0	75 cm	black		NY60DS16
			1	5/0	45 cm	black		NY50DS16-45
			1	5/0	75 cm	black		NY50DS16
			1.5	4/0	45 cm	blue		NY40DS16-45
			1.5	4/0	75 cm	blue		NY40DS16
			2	3/0	75 cm	blue		NY30DS16




Reverse cutting needle – 3/8 circle


NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS19	19 mm		1	5/0	75 cm	black		NY50DS19
			1.5	4/0	45 cm	blue		NY40DS19-45
			1.5	4/0	75 cm	blue		NY40DS19
			2	3/0	75 cm	blue		NY30DS19
			3	2/0	75 cm	blue		NY20DS19
DS24	24 mm		1.5	4/0	75 cm	blue		NY40DS24
			1.5	4/0	45 cm	blue		NY40DS24-45
			2	3/0	75 cm	blue		NY30DS24
			3	2/0	75 cm	blue		NY20DS24
			3.5	0	75 cm	blue		NY00DS24
DS26	26 mm		2	3/0	75 cm	blue		NY30DS26
			3	2/0	75 cm	blue		NY20DS26
DS30	30 mm		2	3/0	75 cm	blue		NY30DS30
			3	2/0	75 cm	blue		NY20DS30
			3.5	0	75 cm	blue		NY00DS30
			4	1	75 cm	blue		NY01DS30
DS36	36 mm		3	2/0	75 cm	blue		NY20DS36
			3.5	0	75 cm	blue		NY00DS36
			4	1	75 cm	blue		NY01DS36
DS40	40 mm		3	2/0	75 cm	blue		NY20DS40
			3.5	0	75 cm	blue		NY00DS40
			4	1	75 cm	blue		NY01DS40
			5	2	90 cm	blue		NY02DS40
DS90	90 mm		3	2/0	100 cm	blue		NY20DS90
			4	1	90 cm	blue		NY01DS90









Reverse cutting needle – 1/2 circle


NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS16	16 mm		1	5/0	75 cm	black		NY50HS16
HS17	17 mm		1.5	4/0	75 cm	blue		NY40HS17
HS36	36 mm		4	1	75 cm	blue		NY01HS36
HS40	40 mm		4	1	75 cm	blue		NY01HS40

Reverse cutting needle – 1/2 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HS50 50 mm 	5	2	75 cm	blue		NY02HS50	

Side cutting spatula – 1/2 circle – double needle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HSP062 6.2 mm 150 µm 	0.2	10/0	30 cm	black		NY10HSP062-D	

Side cutting spatula – 3/8 circle – double needle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DSP060 6 mm 200 µm 	0.2	10/0	30 cm	black		NY10DSP060-D	
DSP062 6.2 mm 150 µm 	0.2	10/0	30 cm	black		NY10DSP062-D	
	0.3	9/0	30 cm	black		NY90DSP062-D	
	0.4	8/0	30 cm	black		NY80DSP062-D	
DSP062 6.2 mm 200 µm 	0.5	7/0	30 cm	black		NY70DSP062-D	
DSP064 6.4 mm 200 µm 	0.2	10/0	30 cm	black		NY10DSP064-D	
DSP080 8 mm 300 µm 	0.7	6/0	45 cm	black		NY60DSP080-D	

Side cutting spatula – 3/8 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DSP062 6.2 mm 150 µm 	0.2	10/0	15 cm	black		NY10DSP062	

Round bodied micro needle – taper point – 3/8 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DSP064 6.4 mm 150 µm 	0.4	8/0	15 cm	black		NY80M-DR064-150	

Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	100	blue		NY40100
2	3/0	100	blue		NY30100
3	2/0	100	blue		NY20100
3.5	0	100	blue		NY00100
4	1	100	blue		NY01100
5	2	75	blue		NY0275
6	3 & 4	50	blue		NY0350