






# LUXAMID



NYLON BLUE



## NONABSORBABLE SURGICAL SUTURE








TYPE	Monofilament
COMPOSITION	Extrusion of polyamide 6.0 or 6.6
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable, gradually encapsulated by connective tissue. The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).
SUTURE COLOUR	Blue, black
RANGE	USP 10/0 - 2
INDICATIONS	Luxamid sutures are intended for use in general soft tissue closing and/or ligation; especially in Luxamid blue: skin closure, general surgery, plastic surgery, gastrointestinal surgery, gynaecology, obstetrics, orthopaedics Luxamid black: plastic surgery, ophthalmic surgery





Round bodied needle – taper point – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16	16 mm		1.5	4/0	75 cm	blue		NY40HR16
HR20	20 mm		2	3/0	75 cm	blue		NY30HR20
HR26	26 mm		3	2/0	75 cm	blue		NY20HR26
HR30	30 mm		2	3/0	75 cm	blue		NY30HR30
			3	2/0	75 cm	blue		NY20HR30
HR40	40 mm		3.5	0	75 cm	blue		NY00HR40
			4	1	75 cm	blue		NY01HR40
			5	2	75 cm	blue		NY02HR40

Round bodied needle – taper point – 1/2 circle – loop								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR40	40 mm		3.5	0	150 cm	blue	loop	NY00HR40-L
			4	1	150 cm	blue	loop	NY01HR40-L
HR50	50 mm		4	1	150 cm	blue	loop	NY01HR50-L

Reverse cutting needle – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12	12 mm		0.5	7/0	75 cm	black		NY70DS12
			0.7	6/0	75 cm	black		NY60DS12
			1	5/0	75 cm	black		NY50DS12
			1.5	4/0	75 cm	blue		NY40DS12
DS16	16 mm		0.7	6/0	45 cm	black		NY60DS16-45
			0.7	6/0	75 cm	black		NY60DS16
			1	5/0	45 cm	black		NY50DS16-45
			1	5/0	75 cm	black		NY50DS16
			1.5	4/0	45 cm	blue		NY40DS16-45
			1.5	4/0	75 cm	blue		NY40DS16
			2	3/0	75 cm	blue		NY30DS16

Reverse cutting needle – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS19	19 mm		1	5/0	75 cm	black		NY50DS19
			1.5	4/0	45 cm	blue		NY40DS19-45
			1.5	4/0	75 cm	blue		NY40DS19
			2	3/0	75 cm	blue		NY30DS19
			3	2/0	75 cm	blue		NY20DS19
DS24	24 mm		1.5	4/0	45 cm	blue		NY40DS24-45
			1.5	4/0	75 cm	blue		NY40DS24
			2	3/0	45 cm	blue		NY30DS24-45
			2	3/0	75 cm	blue		NY30DS24
			3	2/0	75 cm	blue		NY20DS24
			3.5	0	75 cm	blue		NY00DS24
DS26	26 mm		2	3/0	75 cm	blue		NY30DS26
			3	2/0	75 cm	blue		NY20DS26
DS30	30 mm		2	3/0	75 cm	blue		NY30DS30
			3	2/0	75 cm	blue		NY20DS30
			3.5	0	75 cm	blue		NY00DS30
			4	1	75 cm	blue		NY01DS30
DS36	36 mm		2	3/0	75 cm	blue		NY30DS36
			3	2/0	75 cm	blue		NY20DS36
			3.5	0	75 cm	blue		NY00DS36
			4	1	75 cm	blue		NY01DS36
DS40	40 mm		3	2/0	75 cm	blue		NY20DS40
			3	2/0	90 cm	blue		NY20DS40-90
			3.5	0	75 cm	blue		NY00DS40
			4	1	75 cm	blue		NY01DS40
			5	2	90 cm	blue		NY02DS40
DS90	90 mm		3	2/0	100 cm	blue		NY20DS90
			4	1	90 cm	blue		NY01DS90

Reverse cutting needle – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS16	16 mm		1	5/0	75 cm	black		NY50HS16
HS17	17 mm		1.5	4/0	75 cm	blue		NY40HS17