




NONABSORBABLE SURGICAL SUTURE



SUPRAMID



POLYAMID BLACK




TYPE	Twisted and coated multifilament, sizes 5/0 and 6/0 monofilament
COMPOSITION	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable, gradually encapsulated by connective tissue
SUTURE COLOUR	Black
RANGE	USP 6/0 - 2/0
INDICATIONS	Supramid sutures are intended for use in general soft tissue closing and/or ligation; especially in skin closure












 Round bodied needle – taper point – 1/2 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HR16 16 mm 	1	5/0	75 cm	black		SB50HR16	
	1.5	4/0	75 cm	black		SB40HR16	
HR20 20 mm 	1.5	4/0	75 cm	black		SB40HR20	

 Round bodied needle – taper point – straight							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
GR22 22 mm 	1.5	4/0	50 cm	black		SB40GR22	

 Round bodied needle – tapercut point – 3/8 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DRT18 18 mm 	1.5	4/0	75 cm	black		SB40DRT18	
	2	3/0	75 cm	black		SB30DRT18	

 Round bodied needle – tapercut point – 1/2 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HRT16 16 mm 	1	5/0	75 cm	black		SB50HRT16	
HRT18 18 mm 	1.5	4/0	75 cm	black		SB40HRT18	
	2	3/0	75 cm	black		SB30HRT18	

Reverse cutting needle – 3/8 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DS12 12 mm 	0.7	6/0	75 cm	black		SB60DS12	
DS16 16 mm 	0.7	6/0	75 cm	black		SB60DS16	
	1	5/0	45 cm	black		SB50DS16-45	
	1	5/0	75 cm	black		SB50DS16	
	1.5	4/0	75 cm	black		SB40DS16-75	
DS16 16 mm PREMIUM NEEDLE 	1	5/0	75 cm	black		SB50DS16-P	
DS19 19 mm 	1	5/0	75 cm	black		SB50DS19	
	1.5	4/0	45 cm	black		SB40DS19-45	
	1.5	4/0	75 cm	black		SB40DS19	
	2	3/0	75 cm	black		SB30DS19	
DS24 24 mm 	3	2/0	75 cm	black		SB20DS24	

Reverse cutting needle – 1/2 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HS16 16 mm 	1	5/0	75 cm	black		SB50HS16	
HS17 17 mm 	1.5	4/0	75 cm	black		SB40HS17	
HS20 20 mm 	1	5/0	75 cm	black		SB50HS20	
	1.5	4/0	75 cm	black		SB40HS20	
	2	3/0	75 cm	black		SB30HS20	
HS22 22 mm 	1.5	4/0	75 cm	black		SB40HS22	