







Without needles

| NEEDLE | EP  | USP | LENGTH | SUTURE COLOUR | OTHER INFORMATION | ORDER NUMBER |
|--------|-----|-----|--------|---------------|-------------------|--------------|
|        | 3.5 | 0   | 75 cm  | green         | N18               | PE00WN150    |
|        | 4   | 1   | 75 cm  | green         | N18               | PE01WN150    |
|        | 5   | 2   | 150 cm | green         | N18               | PE02WN150    |
|        | 2   | 3/0 | 250 cm | green         | N18               | PE30WN250    |
|        | 3   | 2/0 | 250 cm | green         | N18               | PE20WN250    |
|        |     |     |        |               |                   |              |




Round bodied needle – taper point – 1/2 circle

| NEEDLE   | EP  | USP | LENGTH | SUTURE COLOUR | OTHER INFORMATION | ORDER NUMBER |
|--|-----|-----|--------|---------------|-------------------|--------------|
| HR16 16 mm    | 1.5 | 4/0 | 75 cm  | green         | N18               | PE40HR16     |
| HR20 20 mm    | 2   | 3/0 | 75 cm  | green         | N17               | PE30HR20     |
| HR26 26 mm   | 2   | 3/0 | 75 cm  | green         | N18               | PE30HR26     |
|  | 3   | 2/0 | 75 cm  | green         | N17               | PE20HR26     |
| HR30 30 mm  | 3   | 2/0 | 75 cm  | green         | N18               | PE20HR30     |
|  | 3.5 | 0   | 75 cm  | green         |                   | PE00HR30     |
| HR36 36 mm  | 3.5 | 0   | 75 cm  | green         |                   | PE00HR36     |
|  | 4   | 1   | 75 cm  | green         |                   | PE01HR36     |
| HR40 40 mm  | 4   | 1   | 75 cm  | green         |                   | PE01HR40     |
|  | 5   | 2   | 75 cm  | green         |                   | PE02HR40     |



Round bodied needle – tapercut point – 1/2 circle

| NEEDLE  | EP  | USP | LENGTH | SUTURE COLOUR | OTHER INFORMATION | ORDER NUMBER |
|---|-----|-----|--------|---------------|-------------------|--------------|
| HRT18 18 mm  | 1.5 | 4/0 | 75 cm  | green         | N17               | PE40HRT18    |
|   | 2   | 3/0 | 75 cm  | green         | N17               | PE30HRT18    |



Reverse cutting needle – 3/8 circle

| NEEDLE |       |  | EP  | USP | LENGTH | SUTURE COLOUR | OTHER INFORMATION | ORDER NUMBER |
|--------|-------|--|-----|-----|--------|---------------|-------------------|--------------|
| DS16   | 16 mm |  | 1   | 5/0 | 75 cm  | green         | N20               | PE50DS16     |
|        |       |  | 1.5 | 4/0 | 75 cm  | green         | N20               | PE40DS16     |
| DS19   | 19 mm |  | 1.5 | 4/0 | 75 cm  | green         |                   | PE40DS19     |
|        |       |  | 2   | 3/0 | 75 cm  | green         |                   | PE30DS19     |
| DS24   | 24 mm |  | 2   | 3/0 | 75 cm  | green         |                   | PE30DS24     |
|        |       |  | 3   | 2/0 | 75 cm  | green         |                   | PE20DS24     |
| DS30   | 30 mm |  | 3   | 2/0 | 75 cm  | green         |                   | PE20DS30     |
|        |       |  | 3.5 | 0   | 75 cm  | green         |                   | PE00DS30     |
|        |       |  | 4   | 1   | 75 cm  | green         |                   | PE01DS30     |



Side cutting spatula – 1/4 circle – double needle

| NEEDLE |      |        |  | EP | USP | LENGTH | SUTURE COLOUR | OTHER INFORMATION | ORDER NUMBER |
|--------|------|--------|--|----|-----|--------|---------------|-------------------|--------------|
| VSP080 | 8 mm | 350 µm |  | 1  | 5/0 | 45 cm  | white         | N20               | PE50VSP080-D |