








Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	3.5	0	150 cm	green		PE00WN150
	4	1	150 cm	green		PE01WN150
	5	2	150 cm	green		PE02WN150
	2	3/0	250 cm	green		PE30WN250
	3	2/0	250 cm	green		PE20WN250




Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16 16 mm 	1.5	4/0	75 cm	green		PE40HR16
HR20 20 mm 	2	3/0	75 cm	green		PE30HR20
HR25 25 mm 	3	2/0	75 cm	green	double armed needle	PE20HR25-D
HR26 26 mm 	2	3/0	75 cm	green		PE30HR26
	3	2/0	75 cm	green		PE20HR26
HR30 30 mm 	3	2/0	75 cm	green		PE20HR30
HR36 36 mm 	3.5	0	75 cm	green		PE00HR36
	4	1	75 cm	green		PE01HR36
HR40 40 mm 	4	1	75 cm	green		PE01HR40
	5	2	75 cm	green		PE02HR40



Round bodied needle – tapercut point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18 18 mm 	1.5	4/0	75 cm	green		PE40HRT18
	2	3/0	75 cm	green		PE30HRT18



### Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS16	16 mm		1	5/0	75 cm	green		PE50DS16
			1.5	4/0	75 cm	green		PE40DS16
DS19	19 mm		1.5	4/0	75 cm	green		PE40DS19
			2	3/0	75 cm	green		PE30DS19
DS24	24 mm		2	3/0	75 cm	green		PE30DS24
			3	2/0	75 cm	green		PE20DS24
DS30	30 mm		3.5	0	75 cm	green		PE00DS30



### Side cutting spatula – 1/4 circle – double needle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
VSP080	8 mm	350 µm		1	5/0	45 cm	white	PE50VSP080-D
VSP080	8 mm	430 µm		1.5	4/0	45 cm	white	PE40VSP080-D

### Sutures in cassettes



EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
1.5	4/0	100 m	green		PEG40100
2	3/0	100 m	green		PEG30100
3	2/0	100 m	green		PEG20100
3.5	0	100 m	green		PEG00100
4	1	100 m	green		PEG01100
5	2	75 m	green		PEG0275
6	3 & 4	50 m	green		PEG0350
7	5	40 m	green		PEG0540
8	6	40 m	green		PEG0640
3	2/0	100 m	white		PEW20100
3,5	0	100 m	white		PEW00100
4	1	100 m	white		PEW01100
5	2	100 m	white		PEW02100
6	3 & 4	80 m	white		PEW0380
7	5	60 m	white		PEW0560
8	6	50 m	white		PEW0650
9	7	40 m	white		PEW0740
10	8	30 m	white		PEW0830