









Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	150 cm	violet		PG30WN150
	3	2/0	150 cm	violet		PG20WN150
	3.5	0	150 cm	violet		PG00WN150
	4	1	150 cm	violet		PG01WN150
	5	2	150 cm	violet		PG02WN150
	2	3/0	250 cm	violet	N18	PG30WN250
	3	2/0	250 cm	violet	N18	PG20WN250
	3.5	0	250 cm	violet	N18	PG00WN250
	4	1	250 cm	violet	N18	PG01WN250









Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR13 13 mm 	0.7	6/0	75 cm	violet		PG60HR13
HR17 17 mm 	1	5/0	75 cm	violet		PG50HR17
	1.5	4/0	75 cm	violet	N17	PG40HR17
HR22 22 mm 	1.5	4/0	75 cm	violet		PG40HR22
	3	2/0	75 cm	violet		PG20HR22
	2	3/0	75 cm	violet	N17	PG30HR22
HR26 26 mm 	1.5	4/0	75 cm	violet		PG40HR26
	2	3/0	75 cm	violet		PG30HR26
	3	2/0	75 cm	violet		PG20HR26
HR30 30 mm 	3	2/0	75 cm	violet		PG20HR30
	3.5	0	75 cm	violet	N18	PG00HR30
HR36 36 mm 	2	3/0	75 cm	violet		PG30HR36
	3	2/0	75 cm	violet		PG20HR36
	3.5	0	90 cm	violet		PG00HR36
HR40 40 mm 	3.5	0	90 cm	violet		PG00HR40
	4	1	75 cm	violet	N18	PG01HR40-75
	4	1	90 cm	violet		PG01HR40
	5	2	90 cm	violet		PG02HR40
HR48 48 mm 	3.5	0	90 cm	violet		PG00HR48
	4	1	90 cm	violet		PG01HR48
	5	2	90 cm	violet		PG02HR48







Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12 12 mm 	0.7	6/0	45 cm	violet		PG60DS12
	1	5/0	45 cm	violet		PG50DS12
DS16 16 mm 	1	5/0	75 cm	violet		PG50DS16
	1.5	4/0	75 cm	violet		PG40DS16
	2	3/0	75 cm	violet		PG30DS16
DS19 19 mm 	1	5/0	75 cm	violet		PG50DS19
	1.5	4/0	75 cm	violet		PG40DS19
	1.5	4/0	75 cm	undyed		PG40DS19-U
	2	3/0	75 cm	violet		PG30DS19
	2	3/0	75 cm	undyed		PG30DS19-U
	3	2/0	75 cm	violet		PG20DS19
DS24 24 mm 	2	3/0	75 cm	violet		PG30DS24
	3	2/0	75 cm	violet		PG20DS24
DS30 30 mm 	3	2/0	75 cm	violet		PG20DS30
	3.5	0	75 cm	violet		PG00DS30
	4	1	75 cm	violet		PG01DS30
DS40 40 mm 	3.5	0	75 cm	violet		PG00DS40
	4	1	75 cm	violet		PG01DS40




Reverse cutting needle – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS26 26 mm 	3	2/0	75 cm	violet	N17	PG20HS26
HS36 36 mm 	3.5	0	90 cm	violet	N18	PG00HS36
HS40 40 mm 	3.5	0	90 cm	violet		PG00HS40
	4	1	90 cm	violet		PG01HS40
	5	2	90 cm	violet		PG02HS40
HS48 48 mm 	3.5	0	90 cm	violet		PG00HS48
	4	1	90 cm	violet		PG01HS48
	5	2	90 cm	violet		PG02HS48



Reverse cutting needle – straight

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GS60 60 mm 	3	2/0	75 cm	violet		PG20GS60