


























 Round bodied needle – taper point – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR08	8 mm		0.4	8/0	60 cm	blue	double armed needle	PP80HR08-D
			0.5	7/0	60 cm	blue	double armed needle	PP70HR08-D
HR10	10 mm		0.5	7/0	45 cm	blue		PP70HR10
			0.5	7/0	60 cm	blue	double armed needle	PP70HR10-D
			0.7	6/0	60 cm	blue	double armed needle	PP60HR10-D
HR13	13 mm		0.7	6/0	45 cm	blue		PP60HR13
			0.7	6/0	60 cm	blue	double armed needle	PP60HR13-D
			1	5/0	60 cm	blue	double armed needle	PP50HR13-D
HR15	15 mm		1	5/0	75 cm	blue		PP50HR15
HR16	16 mm		1	5/0	90 cm	blue	double armed needle	PP50HR16-D
			1.5	4/0	75 cm	blue		PP40HR16
			1.5	4/0	90 cm	blue	double armed needle	PP40HR16-D
HR20	20 mm		1.5	4/0	75 cm	blue		PP40HR20
			1.5	4/0	90 cm	blue	double armed needle	PP40HR20-D
			2	3/0	75 cm	blue		PP30HR20
HR25	25 mm		1.5	4/0	75 cm	blue		PP40HR25
			1.5	4/0	90 cm	blue	double armed needle	PP40HR25-D
			2	3/0	75 cm	blue		PP30HR25
			2	3/0	90 cm	blue	double armed needle	PP30HR25-D
			3	2/0	75 cm	blue		PP20HR25
			3	2/0	90 cm	blue	double armed needle	PP20HR25-D
			3.5	0	75 cm	blue		PP00HR25
HR30	30 mm		2	3/0	75 cm	blue		PP30HR30
			3	2/0	75 cm	blue		PP20HR30
			3.5	0	75 cm	blue		PP00HR30
HR35	35 mm		3	2/0	75 cm	blue		PP20HR35
			3.5	0	75 cm	blue		PP00HR35
			4	1	75 cm	blue		PP01HR35
HR40	40 mm		3.5	0	75 cm	blue		PP00HR40
			4	1	90 cm	blue		PP01HR40
			5	2	75 cm	blue		PP02HR40


 Round bodied needle – taper point – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR13	13 mm		0.7	6/0	75 cm	blue		PP60DR13


 Round bodied needle – taper point – 3/8 circle – double armed needle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DR08 8 mm 	0.4	8/0	45 cm	blue		PP80DR08-D
	0.5	7/0	60 cm	blue		PP70DR08-D
DR10 10 mm 	0.5	7/0	60 cm	blue		PP70DR10-D
	0.7	6/0	60 cm	blue		PP60DR10-D
DR13 13 mm 	0.5	7/0	60 cm	blue		PP70DR13-D
	0.7	6/0	60 cm	blue		PP60DR13-D
	1	5/0	75 cm	blue		PP50DR13-D

 Reverse cutting needle – 3/8 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS12 12 mm 	0.7	6/0	45 cm	blue		PP60DS12
	0.7	6/0	75 cm	blue		PP60DS12-75
	1	5/0	45 cm	blue		PP50DS12
DS16 16 mm 	0.7	6/0	45 cm	blue		PP60DS16
	1	5/0	45 cm	blue		PP50DS16
	1	5/0	75 cm	blue		PP50DS16-75
	1.5	4/0	45 cm	blue		PP40DS16-45
	1.5	4/0	75 cm	blue		PP40DS16
	2	3/0	75 cm	blue		PP30DS16
DS19 19 mm 	1	5/0	45 cm	blue		PP50DS19-45
	1	5/0	75 cm	blue		PP50DS19
	1.5	4/0	45 cm	blue		PP40DS19-45
	1.5	4/0	75 cm	blue		PP40DS19
	2	3/0	75 cm	blue		PP30DS19
DS24 24 mm 	1.5	4/0	75 cm	blue		PP40DS24
	2	3/0	75 cm	blue		PP30DS24
	3	2/0	75 cm	blue		PP20DS24
DS30 30 mm 	2	3/0	75 cm	blue		PP30DS30
	3	2/0	75 cm	blue		PP20DS30
	3.5	0	75 cm	blue		PP00DS30
DS36 36 mm 	3.5	2	75 cm	blue		PP00DS36
	4	1	75 cm	blue		PP01DS36
DS40 40 mm 	3	2/0	75 cm	blue		PP20DS40
	3.5	0	75 cm	blue		PP00DS40
	4	1	90 cm	blue		PP01DS40

Reverse cutting needle – straight						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GS60 60 mm 	2	3/0	75 cm	blue		PP30GS60
	3	2/0	75 cm	blue		PP20GS60

Round bodied needle – tapercut point – 1/2 circle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT25 25 mm 	3	2/0	90 cm	blue	double armed needle	PP20HRT25-D
HRT30 30 mm 	2	3/0	90 cm	blue	double armed needle	PP30HRT30-D

Side cutting spatula – 3/8 circle – double needle						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DSP064 6.4 mm 200 µm 	0.2	10/0	30 cm	blue		PP10DSP064-D

Side cutting spatula – straight						
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GSP16 16.15 mm 150 µm 	0.2	10/0	20 cm	blue		PP10GSP16-D