











NONABSORBABLE SURGICAL SUTURE








LUXYLENE



POLYPROPYLENE MONOFILAMENT





TYPE	Monofilament
COMPOSITION	Polypropylene, a polymer of propylene
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable
SUTURE COLOUR	Blue
RANGE	USP 10/0 - 2
INDICATIONS	Luxylene sutures are intended for use in general soft tissue closing and/or ligation; especially in ophthalmic surgery, plastic surgery, skin closure, orthopaedics, hernias (also as a mesh), gastrointestinal surgery, general surgery, gynaecology, obstetrics











 Round bodied needle – taper point – 1/2 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HR10 10 mm 	0.5	7/0	45 cm	blue		PP70HR10	
HR13 13 mm 	0.7	6/0	75 cm	blue		PP60HR13	
HR15 15 mm 	1	5/0	75 cm	blue		PP50HR15	
HR16 16 mm 	1.5	4/0	75 cm	blue		PP40HR16	
HR20 20 mm 	1.5	4/0	75 cm	blue		PP40HR20	
	2	3/0	75 cm	blue		PP30HR20	
HR25 25 mm 	1.5	4/0	75 cm	blue		PP40HR25	
	2	3/0	75 cm	blue		PP30HR25	
	3	2/0	75 cm	blue		PP20HR25	
	3.5	0	75 cm	blue		PP00HR25	
	4	1	75 cm	blue		PP01HR25	
HR30 30 mm 	2	3/0	75 cm	blue		PP30HR30	
	3	2/0	75 cm	blue		PP20HR30	
	3.5	0	75 cm	blue		PP00HR30	
HR35 35 mm 	3	2/0	75 cm	blue		PP20HR35	
	3.5	0	75 cm	blue		PP00HR35	
	4	1	75 cm	blue		PP01HR35	
HR40 40 mm 	3.5	0	75 cm	blue		PP00HR40	
	4	1	90 cm	blue		PP01HR40	
	5	2	75 cm	blue		PP02HR40	


 Round bodied needle – taper point – 1/2 circle – double armed needle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HR08 8 mm 	0.4	8/0	60 cm	blue		PP80HR08-D	
	0.5	7/0	60 cm	blue		PP70HR08-D	
HR10 10 mm 	0.5	7/0	60 cm	blue		PP70HR10-D	
	0.7	6/0	60 cm	blue		PP60HR10-D	
HR13 13 mm 	0.7	6/0	60 cm	blue		PP60HR13-D	
	1	5/0	75 cm	blue		PP50HR13-D	
HR16 16 mm 	1	5/0	90 cm	blue		PP50HR16-D	
	1.5	4/0	90 cm	blue		PP40HR16-D	
HR20 20 mm 	1.5	4/0	90 cm	blue		PP40HR20-D	
HR25 25 mm 	1.5	4/0	90 cm	blue		PP40HR25-D	
	2	3/0	90 cm	blue		PP30HR25-D	
	3	2/0	90 cm	blue		PP20HR25-D	

 Round bodied needle – taper point – 3/8 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DR13 13 mm 	0.7	6/0	75 cm	blue		PP60DR13	

 Round bodied needle – taper point – 3/8 circle – double armed needle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DR08 8 mm 	0.4	8/0	45 cm	blue		PP80DR08-D	
	0.5	7/0	60 cm	blue		PP70DR08-D	
DR10 10 mm 	0.5	7/0	60 cm	blue		PP70DR10-D	
	0.7	6/0	60 cm	blue		PP60DR10-D	
DR13 13 mm 	0.7	6/0	60 cm	blue		PP60DR13-D	
	1	5/0	60 cm	blue		PP50DR13-D	

 Reverse cutting needle – 3/8 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DS12 12 mm 	0.7	6/0	45 cm	blue		PP60DS12	
	0.7	6/0	75 cm	blue		PP60DS12-75	
	1	5/0	45 cm	blue		PP50DS12	
DS16 16 mm 	0.7	6/0	45 cm	blue		PP60DS16	
	1	5/0	45 cm	blue		PP50DS16	
	1	5/0	75 cm	blue		PP50DS16-75	
	1.5	4/0	45 cm	blue		PP40DS16-45	
	1.5	4/0	75 cm	blue		PP40DS16	
	2	3/0	75 cm	blue		PP30DS16	
DS19 19 mm 	1	5/0	45 cm	blue		PP50DS19-45	
	1	5/0	75 cm	blue		PP50DS19	
	1.5	4/0	45 cm	blue		PP40DS19-45	
	1.5	4/0	75 cm	blue		PP40DS19	
	2	3/0	75 cm	blue		PP30DS19	
	3	2/0	75 cm	blue		PP20DS19	
DS24 24 mm 	1.5	4/0	75 cm	blue		PP40DS24	
	2	3/0	75 cm	blue		PP30DS24	
	3	2/0	75 cm	blue		PP20DS24	
DS30 30 mm 	2	3/0	75 cm	blue		PP30DS30	
	3	2/0	75 cm	blue		PP20DS30	
	3.5	0	75 cm	blue		PP00DS30	

Reverse cutting needle – 3/8 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
DS36 36 mm 	3.5	2	75 cm	blue		PP00DS36	
	4	1	75 cm	blue		PP01DS36	
DS40 40 mm 	3	2/0	75 cm	blue		PP20DS40	
	3.5	0	75 cm	blue		PP00DS40	
	4	1	90 cm	blue		PP01DS40	

Reverse cutting needle – straight							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
GS60 60 mm 	2	3/0	75 cm	blue		PP30GS60	
	3	2/0	75 cm	blue		PP20GS60	