

SILK BRAIDED

BLACK

NONABSORBABLE SURGICAL SUTURE


TYPE	Braided multifilament
COMPOSITION	Braided fibres from the cocoon of the silkworm
COATING	Silicone
TISSUE REACTION	Moderate
ABSORPTION	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues
SUTURE COLOUR	Black
RANGE	USP 8/0 - 3&4
INDICATIONS	Silk sutures are intended for use in general soft tissue closing and/or ligation; especially in general surgery, skin closure, ophthalmic surgery, oral surgery, gastrointestinal surgery, gynaecology, obstetrics








Without needles





NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	150 cm	black		SI30WN150
	3	2/0	150 cm	black		SI20WN150
	3.5	0	150 cm	black		SI00WN150
	4	1	150 cm	black		SI01WN150
	2	3/0	13 x 60 cm	black		SI30WN60-13
	3	2/0	13 x 60 cm	black		SI20WN60-13
	3.5	0	13 x 60 cm	black		SI00WN60-13
	4	1	13 x 60 cm	black		SI01WN60-13
	3	2/0	12 x 75 cm	black		SI20WN75-12
	3.5	0	12 x 75 cm	black		SI00WN75-12
	4	1	12 x 75 cm	black		SI01WN75-12
	5	2	12 x 75 cm	black		SI02WN75-12



Without needles – ligapack



NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	3	2/0	250 cm	black		SI20WN250
	5	0	250 cm	black		SI02WN250



 **Round bodied needle – taper point – 1/2 circle**

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR13 13 mm 	1	5/0	45 cm	black		SI50HR13
HR16 16 mm 	1	5/0	75 cm	black		SI50HR16
	1.5	4/0	75 cm	black		SI40HR16
HR20 20 mm 	1	5/0	75 cm	black		SI50HR20
	1.5	4/0	75 cm	black		SI40HR20
	2	3/0	75 cm	black		SI30HR20
	3	2/0	75 cm	black		SI20HR20
HR26 26 mm 	1.5	4/0	75 cm	black		SI40HR26
	2	3/0	75 cm	black		SI30HR26
	3	2/0	75 cm	black		SI20HR26
	3.5	0	75 cm	black		SI00HR26
	4	1	75 cm	black		SI01HR26

 Round bodied needle – taper point – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR30	30 mm		2	3/0	75 cm	black		SI30HR30
			3	2/0	75 cm	black		SI20HR30
			3.5	0	75 cm	black		SI00HR30
			4	1	75 cm	black		SI01HR30
HR37	37 mm		3.5	0	75 cm	black		SI00HR37
			4	1	75 cm	black		SI01HR37
			5	2	75 cm	black		SI02HR37
HR40	40 mm		3.5	0	75 cm	black		SI00HR40
			4	1	75 cm	black		SI01HR40
			5	2	75 cm	black		SI02HR40









 Round bodied needle – taper point – straight								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
GR22	22 mm		1.5	4/0	50 cm	black		SI40GR22





 Round bodied needle – tapercut point – 3/8 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DRT18	18 mm		1.5	4/0	75 cm	black		SI40DRT18
			2	3/0	75 cm	black		SI30DRT18

 Round bodied needle – tapercut point – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18	18 mm		1.5	4/0	75 cm	black		SI40HRT18
			2	3/0	75 cm	black		SI30HRT18



Reverse cutting needle – 3/8 circle

NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS8	8 mm		0.7	6/0	45 cm	black		SI60DS08
DS12	12 mm		0.7	6/0	45 cm	black		SI60DS12
			1	5/0	45 cm	black		SI50DS12
DS16	16 mm		0.7	6/0	75 cm	black		SI60DS16
			1	5/0	75 cm	black		SI50DS16
			1.5	4/0	75 cm	black		SI40DS16
			2	3/0	75 cm	black		SI30DS16
DS19	19 mm		1	5/0	75 cm	black		SI50DS19
			1.5	4/0	45 cm	black		SI40DS19-45
			1.5	4/0	75 cm	black		SI40DS19
			2	3/0	75 cm	black		SI30DS19
			3	2/0	75 cm	black		SI20DS19
DS24	24 mm		2	3/0	75 cm	black		SI30DS24
			3	2/0	75 cm	black		SI20DS24
			3.5	0	75 cm	black		SI00DS24
DS30	30 mm		1.5	4/0	75 cm	black		SI40DS30
			2	3/0	75 cm	black		SI30DS30
			3	2/0	75 cm	black		SI20DS30
			3.5	0	75 cm	black		SI00DS30
			4	1	75 cm	black		SI01DS30
DS37	37 mm		2	3/0	75 cm	black		SI30DS37
			3	2/0	75 cm	black		SI20DS37
			3.5	0	75 cm	black		SI00DS37
			4	1	75 cm	black		SI01DS37
			5	2	75 cm	black		SI02DS37
DS45	45 mm		3.5	0	75 cm	black		SI00DS45
			4	1	75 cm	black		SI01DS45
			5	2	75 cm	black		SI02DS45

Reverse cutting needle – 1/2 circle								
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS17	17 mm		1.5	4/0	75 cm	black		SI40HS17
			2	3/0	75 cm	black		SI30HS17
HS20	20 mm		1.5	4/0	75 cm	black		SI40HS20
			2	3/0	75 cm	black		SI30HS20
HS22	22 mm		1.5	4/0	75 cm	black		SI40HS22
			2	3/0	75 cm	black		SI30HS22
			3	2/0	75 cm	black		SI20HS22
HS26	26 mm		2	3/0	75 cm	black		SI30HS26
			3	2/0	75 cm	black		SI20HS26