




NONABSORBABLE SURGICAL SUTURE






STEEL



STAINLESS STEEL MONOFILAMENT

TYPE	Monofilament
COMPOSITION	Iron, nickel and chromium alloy
COATING	None
TISSUE REACTION	Minimal
ABSORPTION	Non absorbable
SUTURE COLOUR	Metal
RANGE	USP 1 - 7
INDICATIONS	Steel sutures are intended for use in closure of sternum, orthopaedics, abdominal wall closure and hernia repair



 Round bodied needle – tapercut point – 1/2 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HRT40 40 mm 	6	3 & 4	4 x 45 cm	metal		ST03HRT40	
HRT48 48 mm 	7	5	4 x 45 cm	metal		ST05HRT48	
	9	7	4 x 45 cm	metal		ST07HRT48	

 Reverse cutting needle – 1/2 circle							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
HS40 40 mm 	4	1	4 x 45 cm	metal		ST01HS40	
	5	2	4 x 45 cm	metal		ST02HS40	
HS48 48 mm 	6	3 & 4	4 x 45 cm	metal		ST03HS48	
HS50 50 mm 	7	5	4 x 45 cm	metal		ST05HS50	
	8	6	4 x 45 cm	metal		ST06HS50	
HS58 58 mm 	7	5	4 x 45 cm	metal		ST05HS58	
	8	6	4 x 45 cm	metal		ST06HS58	

 Reverse cutting needle – 1/2 circle – PATELLA SET							
NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER	
NH120 120 mm 	9	7	60 cm	metal	PATELLA SET	ST07HS120	